

TLG® PLASTIC LABORATORY BOTTLES (HDPE)



TLG® Plastic Laboratory Bottles from Sati International meet the requirements of the most demanding applications. These plastic bottles are highly durable and strong for a safe, shatterproof alternative to glass bottles for collection, storage, or shipping of samples, media, and reagents.

Narrow Mouth High-Density Polyethylene (HDPE)

- Safe for freezer storage down to -100°C
- Resistant to most corrosives

Catalog Number	Material	Capacity	Cap Size
CS-C1615030	HDPE	30 ML	20 mm
CS-C1615060	HDPE	60 ML	20 mm
CS-C1615125	HDPE	125 ML	24 mm
CS-C1615250	HDPE	250 ML	24 mm
CS-C1615500	HDPE	500 ML	24 mm
CS-C161501L	HDPE	1000 ML	38 mm

Wide Mouth High-Density Polyethylene (HDPE)

- Safe for freezer storage down to -100°C
- Resistant to most corrosives

Catalog Number	Material	Capacity	Cap Size
CS-C1621030	HDPE	30 ML	28 mm
CS-C1621060	HDPE	60 ML	28 mm
CS-C1621125	HDPE	125 ML	38 mm
CS-C1621250	HDPE	250 ML	43 mm
CS-C1621500	HDPE	500 ML	43 mm
CS-C16211L	HDPE	1000 ML	63 mm
CS-C16212L	HDPE	2000 ML	83 mm

TLG® PLASTIC LABORATORY BOTTLES (LDPE)



TLG® Plastic Laboratory Bottles from Sati International meet the requirements of the most demanding applications. These plastic bottles are highly durable and strong for a safe, shatterproof alternative to glass bottles for collection, storage, or shipping of samples, media, and reagents.

Narrow Mouth Low-Density Polyethylene (LDPE)

- More clear than HDPE
- Flexible with excellent impact resistance

Catalog Number	Material	Capacity	Cap Size
CS-C1614030	LDPE	30 ML	20 mm
CS-C1614060	LDPE	60 ML	20 mm
CS-C1614125	LDPE	125 ML	24 mm
CS-C1614250	LDPE	250 ML	24 mm
CS-C1614500	LDPE	500 ML	28 mm
CS-C161401L	LDPE	1000 ML	38 mm

Wide Mouth Low-Density Polyethylene (LDPE)

- More clear than HDPE
- Flexible with excellent impact resistance

Catalog Number	Material	Capacity	Cap Size
CS-C1623030	LDPE	30 ML	28 mm
CS-C1623060	LDPE	60 ML	28 mm
CS-C1623125	LDPE	125 ML	38 mm
CS-C1623250	LDPE	250 ML	43 mm
CS-C1623500	LDPE	500 ML	53 mm
CS-C16231L	LDPE	1000 ML	63 mm
CS-C16232L	LDPE	2000 ML	83 mm

TLG® PLASTIC LABORATORY BOTTLES (PP)



TLG® Plastic Laboratory Bottles from Sati International meet the requirements of the most demanding applications. These plastic bottles are highly durable and strong for a safe, shatterproof alternative to glass bottles for collection, storage, or shipping of samples, media, and reagents.

Narrow Mouth, Polypropylene (PP)

- Suitable for autoclaving at 121°C
- Improved chemical resistance and clarity
- Resistant to most corrosives and suitable for a wide variety of laboratory uses
- Approved for use from -100 °C to 121°C

Catalog Number	Material	Capacity	Cap Size
CS-C1616030	PP	30 ML	20 mm
CS-C1616060	PP	60 ML	20 mm
CS-C1616125	PP	125 ML	24 mm
CS-C1616250	PP	250 ML	24 mm
CS-C1616500	PP	500 ML	24 mm
CS-C161601L	PP	1000 ML	38 mm

Wide Mouth, Polypropylene (PP)

- Suitable for autoclaving at 121°C
- Improved chemical resistance and clarity
- Resistant to most corrosives and suitable for a wide variety of laboratory uses
- Approved for use from -100 °C to 121°C

Catalog Number	Material	Capacity	Cap Size
CS-C1622030	PP	30 ML	28 mm
CS-C1622060	PP	60 ML	28 mm
CS-C1622125	PP	125 ML	38 mm
CS-C1622250	PP	250 ML	43 mm
CS-C1622500	PP	500 ML	43 mm
CS-C16221L	PP	1000 ML	63 mm