

# LABWARE PLASTIC



 1 (833) 636-7850

 [WWW.SATIINTERNATIONAL.COM](http://WWW.SATIINTERNATIONAL.COM)

**Chem Science Inc**

965 Reverchon St-Laurent (Quebec) Canada H4T 4L2

## Plastic Beakers without Handle

- Manufacture in state of the art molding facility using high purity virgin polypropylene (PP).
- Used for general laboratory use.
- Tapered stacking style saves space and makes it easy to handle.
- All sizes have to read silk screened graduations.
- Precise pouring spouts.
- Can be autoclaved.



Catalog No.	Capacity	Height	Pack Type
CS-P050610	10 ML	32 mm	12
CS-P050609	25 ML	42 mm	12
CS-P050601	50 ML	51 mm	12
CS-P050602	100 ML	62 mm	12
CS-P050601A	150 ML	76.2 mm	12
CS-P050603	250 ML	82 mm	12
CS-P050612	400 ML	82 mm	12
CS-P050604	500 ML	100 mm	12
CS-P050611	600 ML	100 mm	12
CS-P050605	1000 ML	127 mm	12
CS-P050606	2000 ML	158 mm	4
CS-P050607	5000 ML	263 mm	2

## Plastic Beakers with Handle

- Manufacture using high purity virgin polypropylene (PP).
- Used for general laboratory use.
- Rectangular handle for easy carrying and pouring.
- Graduation color - BLUE print
- Precise pouring spouts.
- Can be autoclaved.



Catalog No.	Capacity	Height	Pack Type
CS-P050803	250 mL	94 mm	24
CS-P050804	500 mL	128 mm	24
CS-P050805	1000 mL	150 mm	24
CS-P050806	2000 mL	192 mm	12
CS-P050809	3000 mL	254 mm	12
CS-P050807	5000 mL	263 mm	12
CS-P050808	10000 mL	330 mm	12

## Measuring Cylinder, PP

- Molded in large, easy to read embossed graduations.
- Molded in highest quality translucent virgin polypropylene.
- Cylinders have generous pour spout for smooth flow of liquid during pouring.
- Hexagonal shape gives stability to the base.
- Can be autoclaved.
- Meet accuracy requirements of class B under "DIN 1261/ISO 6706"



Catalog No.	Capacity	Sub Division	Tolerance	Height	Pack Type
CS-P050301	10 ML	0.2 ML	±0.2 ML	138 mm	12
CS-P050302	25 ML	0.5 ML	±0.5 ML	170 mm	12
CS-P050303	50 ML	1.0 ML	±1.0 ML	198 mm	12
CS-P050304	100 ML	1.0 ML	±1.0 ML	257 mm	10
CS-P050305	250 ML	2.0 ML	±2.0 ML	314 mm	6
CS-P050306	500 ML	5.0 ML	±5.0 ML	344 mm	4
CS-P050307	1000 ML	10.0 ML	±10.0 ML	412 mm	2
CS-P050308	2000 ML	20.0 ML	±20.0 ML	482 mm	2

## Erlenmeyer Flask, Wide Neck, PP

- Translucent polypropylene flasks are manufactured with blue graduations for reference.
- The flasks have a screw neck and are supplied with a PP cap.



Catalog No.	Capacity	O.D.	Pack Type
CS-C266025	25 ML	52 mm	12
CS-C266050	50 ML	52 mm	12
CS-C266100	100 ML	52 mm	12
CS-C266125	125 ML	52 mm	12
CS-C266250	250 ML	85 mm	12
CS-C266500	500 ML	108 mm	6
CS-C26601L	1000 ML	135 mm	4
CS-C26602L	2000 ML	135 mm	2



# Separatory Funnel, Pear Shape, Screw Cap, PP



Catalog No.	Capacity	Approx O.D. x Height	Closure	Stem Length	Pack Type
CS-C271125	125 ML	216x67 mm	24.7	65 mm	6
CS-C271250	250 ML	245X81 mm	32.5	65 mm	6
CS-C271500	500 ML	304X101 mm	32.5	65 mm	6
CS-C2711L	1000 ML	325X134 mm	24.7	65 mm	6

Supplying Specialist of  
Laboratory Glassware & Consumables to  
**Scientist and Researchers**



## Sati International Scientific Inc.

965 Reverchon St-Laurent ,(Quebec) Canada H4T 4L2  
Toll Free : 1 (833) 636-7850 | [sales@satiinternational.com](mailto:sales@satiinternational.com)  
[www.satiinternational.com](http://www.satiinternational.com)