

PYCNOMETERS GAY LUSSAC, CALIBRATED



About OUR PRODUCT

Manufactured from ASTM 438 type 1 Borosilicate 3.3 Glass. Precisely balanced at 20° Centigrade, while the joints are 10/18 size. The capillary stoppers are not interchangeable since the final adjustment for volume is made on the stopper. Adjusted bottles are calibrated to ASTM Test Reference D-369. Just fill with your sample and weigh. Each bottle and stopper convey the same recognizing number. Comes with Individual Calibration Certificate.

Catalog No.	Capacity	Joint Size	Tolerance	Pack Type
268.202.00	2 mL	10/18	±0.03 ML	Pack of 2
268.202.0A	5 mL	10/18	±0.03 ML	Pack of 2
268.202.01	10 mL	10/18	±0.04 ML	Pack of 2
268.202.02	25 mL	10/18	±0.05 ML	Pack of 2
268.202.03	50 mL	10/18	±0.08 ML	Pack of 2
268.202.04	100 mL	10/18	±0.10 ML	Pack of 2

Sati International Scientific Inc.

965 Reverchon St-Laurent ,(Quebec) Canada H4T 4L2
www.satiinternational.com | sales@satiinternational.com

TOLL FREE:

1 (833) 636-7850