

WIDE MOUTH

SAFETY WASH BOTTLE



FEATURES

- ✓ Wide neck for easy filling
- ✓ Color – coded lids and contents
- ✓ Made up of LDPE with PP closure
- ✓ Dispensing nozzles are carefully molded to deliver a straight stream for precision dispensing and safety
- ✓ CAS numbers on every bottle
- ✓ It display the chemical name and formula for easy identification
- ✓ The Globally Harmonized system labeling is also present on the bottles
- ✓ Signal word to identify the severity of the hazard elements
- ✓ It also displays the list of recommended protective equipments need while working with bottles

Catalog Number	For use with	Capacity	Stopper Color	Pack Size
CS-C460105E	Acetone	500 mL	Red	6
CS-C460115E	Water	500 mL	Blue	6
CS-C460120E	Ethanol	500 mL	Orange	6
CS-C460130E	Isopropanol	500 mL	Yellow	6
CS-C460140E	Methanol	500 mL	Green	6
CS-C460150E	Ethyl Acetate	500 mL	Green	6
CS-C460160E	Dichloromethane	500 mL	Yellow	6
CS-C460180E	Methyl-Ethyl Ketone	500 mL	Green	6
CS-C460190E	Acetonitrile	500 mL	Yellow	6
CS-C460200E	Sodium Hypochlorite	500 mL	Yellow	6