



**LIMITED
TIME
OFFER**

REUSABLE BOROSILICATE GLASS PETRI DISHES

These Glass Petri Dishes come with beaded edges and can withstand repeated sterilization. The beaded edges provide good mechanical strength and also serve the purpose of equally spacing the sidewalls of the bottom and cover, thus reducing the capillary action of condensed moisture on the sides. They are not affected chemically or thermally by any of the methods commonly employed in laboratories where sterilization is a major factor in routine or in specialized work.

Catalog No.	Description	Diameter	Height	Pack Size	Price
CS-PD-6015	Petri Dish, Complete Set, 60x15 mm	60 mm	15 mm	12	\$ 84.90
CS-PD-10020	Petri Dish, Complete Set, 100x20 mm	100 mm	20 mm	12	\$ 66.90
CS-PD-15020	Petri Dish, Complete Set, 150x20mm	150 mm	20 mm	12	\$ 176.45

Sati International Scientific Inc

965 Reverchon St-Laurent (Quebec) Canada H4T 4L2

Website : www.satiinternational.com

Toll Free : 1 (833) 636-7850

Email : sales@satiinternational.com