

POWDER FUNNEL

POLYPROPYLENE

Transfer powders and thick creams with TLG® Polypropylene Powder Funnels into narrow bottles without spillage. Made from high quality virgin, laboratory grade, polypropylene to ensure excellent results. Available in different sizes for different container sizes and volumes.

- Polypropylene is break-resistant and chemically compatible with a wide variety of lab reagents.
- Suitable for general lab use, it is a durable replacement for glass funnels.
- Parallel stem minimises bridging of powder.
- External ribs prevent air lock, ensuring effective powder transfer.
- Autoclavable.



LIMITED TIME
OFFER
TILL STOCK LAST

| Catalog No. | Capacity | Top I.D. | Stem O.D. | Case Size | Promotional Price |
|-------------|----------|----------|-----------|-----------|-------------------|
| CS-C268065 | 50 mL | 65 mm | 18 mm | 36 | Sold Out |
| CS-C268080 | 100 mL | 80 mm | 22 mm | 36 | \$ 90.90/Case |
| CS-C268100 | 225 mL | 100 mm | 26 mm | 24 | Sold Out |
| CS-C268150 | 750 mL | 150 mm | 36 mm | 24 | \$ 98.30/Case |

Toll Free:
1 (833) 636-7850

Email : sales@satiinternational.com

Website : www.satiinternational.com