



LIMITED TIME
OFFER
Till Stock Last

LABORATORY BOTTLES

POLYPROPYLENE (PP)

TLG® Plastic Laboratory Bottles from Sati International meet the requirements of the most demanding applications. These plastic bottles are highly durable and strong for a safe, shatterproof alternative to glass bottles for collection, storage, or shipping of samples, media and reagents.

- ▶ Caps provide the ultimate leakproof protection
- ▶ Provides excellent chemical Resistant for a variety of lab applications
- ▶ Translucence provides good clarity for viewing contents
- ▶ Autoclavable

Clear, Narrow Mouth Polypropylene (PP) Bottles

Catalog No.	Material	Cap Size	Capacity	Case Size	Price/Case
CS-C1616030	PP	20 mm	30 ML	72	\$ 58.90
CS-C1616060	PP	20 mm	60 ML	72	\$ 58.90
CS-C1616125	PP	24 mm	125 ML	72	\$ 75.30
CS-C1616250	PP	24 mm	250 ML	72	\$ 103.90
CS-C1616500	PP	28 mm	500 ML	48	\$ 115.90

Clear, Wide Mouth Polypropylene (PP) Bottles

Catalog No.	Material	Cap Size	Capacity	Case Size	Price/Case
CS-C1622030	PP	28 mm	30 ML	72	\$ 98.90
CS-C1622060	PP	28 mm	60 ML	72	\$ 98.90
CS-C1622125	PP	38 mm	125 ML	72	\$ 105.90
CS-C1622250	PP	43 mm	250 ML	72	\$ 149.90
CS-C1622500	PP	43 mm	500 ML	48	\$ 153.90