



Mortar and Pestle

Available in Glass and Porcelain

These tools are used to grind or crush substances into fine powders or pastes. Mortars are bowl-shaped Containers and pestles are heavy, rod shaped tools. Sati International provides the exclusive quality sold as sets but also available separately.

TLG Porcelain Pestles

- ▶ For hand comfort in grinding
- ▶ Unglazed on grinding surface
- ▶ Autoclavable
- ▶ Oversized

Catalog No.	Height	For use with (Mortar)	Pack Size
529.303.01P	114 mm	50 ml, 65 ml, 145 ml	6
529.303.04P	180 mm	275 mL	4
529.303.05P	194 mm	400 mL	4
529.303.06P	222 mm	750 mL	2
529.303.07P	235 mm	1100 mL	2
529.303.08P	262 mm	1900 mL	2



TLG Porcelain Mortar

- ▶ Color : White
- ▶ Deep Form
- ▶ With Lip
- ▶ Available in Various Capacities

Catalog No.	Capacity	Height	Dia. (Metric) Outer	Pack Size
529.303.01M	50 mL	47 mm	70 mm	6
529.303.02M	65 mL	52 mm	80 mm	6
529.303.03M	145 mL	69 mm	90 mm	6
529.303.04M	275 mL	70 mm	115 mm	4
529.303.05M	400 mL	81 mm	130 mm	4
529.303.06M	750 mL	110 mm	160 mm	2
529.303.07M	1100 mL	115 mm	186 mm	2
529.303.08M	1900 mL	130 mm	210 mm	2



Mortar and Pestle, Glass

- ▶ Made of Clear and heavy glass
- ▶ Mortars are made of heavy glass with heavy molded base and lip for pouring.
- ▶ Flat bottom for stability
- ▶ Offer a clear view of materials

Catalog No.	Capacity	Diameter (mm)	Height	Pack Size
CS-MPG-1	1 Ounce	70 mm	60 mm	12
CS-MPG-2	2 Ounce	80 mm	68 mm	12
CS-MPG-4	4 Ounce	100 mm	93 mm	12
CS-MPG-8	8 Ounce	110 mm	98 mm	12
CS-MPG-16	16 Ounce	135 mm	120 mm	Each

