


STERILE & NON STERILE

Centrifuge

Tubes

With Flat Cap

Available in Bulk and Rack Pack



Product Specification

- ✓ Leak proof
- ✓ G-Force up to 20000
- ✓ Temperature from -80° to 121°C
- ✓ High capacity centrifuge rotors
- ✓ Non cytotoxic
- ✓ DNase-/RNase Free
- ✓ Pyrogen Free
- ✓ Large writing surface on body
- ✓ Clear transparent surface
- ✓ Clear and Visible graduations on tube body

Catalog No.	Description	Capacity	Type	Pack Type
CS-CT-15-B	15 ML Centrifuge Conical Tube With Flat Cap, Bulk	15 ML	Non-Sterile	Pack of 500
CS-CT-15-B-MF	15 ML Metal Free Centrifuge Tube With Flat Cap, Bulk	15 ML	Metal Free	Pack of 500
CS-CT-15-B-S	15 ML Centrifuge Conical Tube With Flat Cap, Sterile, Bulk	15 ML	Sterile	Pack of 500
CS-CT-15-R-S	15 ML Centrifuge Conical Tube With Flat Cap, Rack, Sterile	15 ML	Sterile	Pack of 500
CS-CT-15-R-S-MF	15 ML Metal Free Centrifuge Tube With Flat Cap, Sterile, Bulk	15 ML	Metal Free	Pack of 500
CS-CT-50-B	50 ML Centrifuge Conical Tube With Flat Cap, Bulk	50 ML	Non-Sterile	Pack of 500
CS-CT-50-B-MF	50 ML Metal Free Centrifuge Tube With Flat Cap, Bulk,	50 ML	Metal Free	Pack of 500
CS-CT-50-B-S	50 ML Centrifuge Conical Tube With Flat Cap, Bulk, Sterile	50 ML	Sterile	Pack of 500
CS-CT-50-R-S	50 ML Centrifuge Conical Tube With Flat Cap, Rack, Sterile	50 ML	Sterile	Pack of 500
CS-CT-50-K-B	50 ML Self Standing Centrifuge Tube With Flat Cap, Bulk	50 ML	Non-Sterile	Pack of 500
CS-CT-50-K-B-S	50mL Self Standing Centrifuge Tubes With Flat Cap, Bulk, Sterile	50 ML	Self-Standing, Sterile	Pack of 500