


# SEROLOGICAL DISPOSABLE PIPETTES

**STERILE**



- ▶ Made of crystal-clear polystyrene
- ▶ Compatible with all popular pipettors
- ▶ More volume due to reverse (negative) graduations
- ▶ Choice of working method due to negative scaling
- ▶ Simple volume identification with international color code for differentiation on each individual blister pack
- ▶ Available in separate sterile packaging
- ▶ Simple to open and anti-static

Catalog No.	Description	Capacity	Type	Pack Type
CS-85-7125	1 ml Serological Pipettes, Sterile, Individually wrapped	1 ml	Sterile	1800 Pcs/Case
CS-85-7126	2 ml Serological Pipettes, Sterile, Individually wrapped	2 ml	Sterile	1800 Pcs/Case
CS-85-7127	5 ml Serological Pipettes, Sterile, Individually wrapped	5 ml	Sterile	1200 Pcs/Case
CS-85-7128	10 ml Serological Pipettes, Sterile, Individually wrapped	10 ml	Sterile	1200 Pcs/Case
CS-85-7129	25 ml Serological Pipettes, Sterile, Individually wrapped	25 ml	Sterile	600 Pcs/Case
CS-85-7130	50 ml Serological Pipettes, Sterile, Individually wrapped	50 ml	Sterile	450 Pcs/Case