



PS (POLY STYRENE) PETRI DISHES

These PS Petri Dishes are suitable for the microbiological application.

1. Vented Petri dishes have elevated lid which allows for good gas exchange and encourages evaporation
2. This design is useful when the poured dishes are to be used as soon as possible
3. Gamma sterilization is as per ISO1117 norms
4. The design of the dish maintains sterility because the particles would have to go up and over the dish wall to get inside which is very rare in normal airflow.

| Catalog No. | Description | Size | Case Type |
|-------------|---------------------------------------|-----------|-------------|
| CS-PS-65 | PS Petri dish, size 60x15, sterile | 60x15 mm | Case of 500 |
| CS-PS-95 | PS Petri dish, size 90x15mm , sterile | 90x15 mm | Case of 500 |
| CS-PS-105 | PS Petri dish , size 100x15, sterile | 100x15 mm | Case of 500 |