

TOOLED NECK VIALS



TLG® Tooled-Neck Sample Vials are excellent for packaging drug products. Necks are tooled for an accurate fit with polyethylene closures. Vials and closures are supplied in both pack and case quantities.

- Manufactured from borosilicate glass
- Polyethylene closure
- Two piece closure
- Ensures tight and accurate fit with the closures.

Catalog No.	Description of Vial	Approx O.D.x Height	Pack Size
CS-24-1	Clear Glass Tooled-Neck Vials with Polyethylene Closures, Open Bottom Closure , Capacity: 1 dr. (3.7mL)	15 x 45 mm	144
CS-24-2	Clear Glass Tooled-Neck Vials with Polyethylene Closures, Open Bottom Closure,Capacity: 2 dr. (7.4mL)	19 x 48 mm	144
CS-24-3	Clear Glass Tooled-Neck Vials with Polyethylene Closures, Closed Bottom Closure, Capacity: 3 dr. (11.1mL)	21 x 50 mm	144
CS-24-4	Clear Glass Tooled-Neck Vials with Polyethylene Closures, Closed Bottom Closure, Capacity: 4 dr. (14.8mL)	25 x 52 mm	144
CS-24-5	Clear Glass Tooled-Neck Vials with Polyethylene Closures, Closed Bottom Closure, Capacity: 5 dr. (18.5mL)	27 x 55 mm	144
CS-24-7	Clear Glass Tooled-Neck Vials with Polyethylene Closures, Closed Bottom Closure, Capacity: 7 dr. (25.9mL)	29 x 65 mm	144