

GLASS CENTRIFUGE TUBES



TLG® Disposable Glass Conical Centrifuge Tubes, Ungraduated, without Screw Cap

To make pipetting faster and easier, these disposable glass pipets include a 'To Deliver' calibration with blow-out and long, thin, tapered tips. All pipets include negative graduations and fire-polished tips to ensure burr-free, uniform flow rates. Nonsterile unplugged and in large quantities. The standard is ASTM E-714. Calibrated "to deliver" (TD) with blow-out. Pipetting is faster and easier with TLG® Disposable Bulk Glass Serological Pipets' long, thin, tapered tips.



Catalog No.	Description	Capacity	Case Size
CS-C5-ST	TLG® 5 mL Disposable Glass Conical Centrifuge Tubes, without Screw Cap	5 ml	Case of 125
CS-C10-ST	TLG® 10 mL Disposable Glass Conical Centrifuge Tubes, without Screw Cap	10 ml	Case of 125
CS-C15-ST	TLG® 15 mL Disposable Glass Conical Centrifuge Tubes, without Screw Cap	15 ml	Case of 125
CS-C50-ST	TLG® 50 mL Disposable Glass Conical Centrifuge Tubes, without Screw Cap	50 ml	Case of 125

Round-Bottom Glass Centrifuge Tubes

The round bottom, ungraduated centrifuge tubes have a pourout spout for easy decanting.



Catalog No.	Stopper	Capacity	Length	Diameter	Closure Included	Case Size
CS-558-2450ST	24-410	50 ml	122 mm	29 mm	Yes	Case of 12
CS-558-50GS	13	50 ml	133 mm	29 mm	Yes	Case of 12