

## Metal Wire Racks for Centrifuge Tubes


- ▶ Wire Racks with a Chemically Resistant HDPE Coating
- ▶ Racks that are strong enough for transport or sample handling
- ▶ Fits tubes and bottles with diameters ranging from 16 to 20mm, 25-30mm, 55-60 mm.
- ▶ Temperatures as low as -40°C can be tolerated.
- ▶ To keep tubes secure, it has square wells and three deck levels.
- ▶ A slotted bottom provides additional support.



Catalog No.	Color	Position	Diameter (Metric) Well	Dimensions	Pack Size
CS-HS23224-B	Blue	6	55-60 mm	8.1 x 5.3 x 4.1 in (20.5 x 13.5 x 10.4 cm)	Each
CS-HS232461-B	Blue	8	25-30 mm	5.8 x 3.2 x 3.5 in (14.7 x 8 x 8.9 cm)	Case of 96
CS-HS23215-B	Blue	15	16-20 mm	5 x 3 x 3 in (12.7 x 7.6 x 7.6 cm)	Case of 96
CS-120764-B	Blue	40	16-20 mm	9.5 x 3.9 x 3.1 in (24 x 9 x 8 cm)	Case of 48
CS-120770-B	Blue	48	25-30 mm	10.9 x 8.3 x 3.9 in (27.7 x 20.9 x 9.8 cm)	Each
CS-120765-B	Blue	80	16-20 mm	9.5 x 7.6 x 3.1 in (24 x 19 x 8 cm)	Case of 24

## Metal Wire Rack for Test Tubes

- ▶ Wire Racks with a Chemically Resistant HDPE Coating
- ▶ Racks that are strong enough for transport or sample handling
- ▶ Fits tubes and bottles with diameters ranging from 10 to 40 mm.
- ▶ Temperatures as low as -40°C can be tolerated.
- ▶ To keep tubes secure, it has square wells and three deck levels.
- ▶ A slotted bottom provides additional support.



Catalog No.	Color	Position	Diameter (Metric) Well	Dimensions	Pack Size
CS-120769-B	Blue	24	20-25 mm	10.7 x 4.1 x 3.2 in (27.3 x 10.5 x 8 cm)	Each
CS-120771-B	Blue	24	30-40 mm	14.1 x 5.5 x 3.9 in (35.7 x 13.9 x 9.8 cm)	Case of 12
CS-120760-B	Blue	40	13-16 mm	7.9 x 3.3 x 2.4 in (20 x 8.3 x 5.9 cm)	Each
CS-120766-B	Blue	40	20-24 mm	10.2 x 4.2 x 3.2 in (26 x 10.7 x 8 cm)	Case of 24
CS-120767-B	Blue	40	20-25 mm	11.4 x 4.7 x 3.1 in (29 x 11.8 x 8 cm)	Case of 24
CS-120761-B	Blue	48	13-16 mm	9.4 x 3.3 x 2.4 in (23.7 x 12.7 x 5.9 cm)	Case of 56
CS-120772-B	Blue	48	30-40 mm	14 x 10.6 x 3.8 in (35 x 26 x 9.8 cm)	Case of 8
CS-120089-B	Blue	72	10-13 mm	8.1 x 4.1 x 2.6 in (20.5 x 10.5 x 6.5 cm)	Case of 48
CS-HS23072-B	Blue	72	13-16 mm	9.5 x 5 x 2.6 in (24.1 x 18.8 x 6.5 cm)	Case of 24
CS-120768-B	Blue	80	20-25 mm	11.4 x 9.1 x 3.1 in (29 x 23 x 8 cm)	Case of 12
CS-120759-B	Blue	108	10-13 mm	8 x 6.4 x 2.6 in (20.4 x 16.1 x 6.5 cm)	Case of 36
CS-120763-B	Blue	108	13-16 mm	9.3 x 5 x 2.6 in (23.7 x 12.7 x 6.5 cm)	Case of 24