

# ERLENMEYER FLASK

## NARROW AND WIDE NECK


Sati International provides Erlenmeyer Flasks that can be used for storing, shaking, and titration. They are made up of borosilicate glass, flat bottom design with durable white enamel scale.

Limited Time  
**OFFER**  
Till Stock Last



## Erlenmeyer Flask, Wide Neck With Graduation

- Made of borosilicate glass 3.3.
- Feature durable white ceramic enamel scale
- Which indicate approximate volume for measuring and mixing.
- Reinforced rim for easy pouring and also avoids any breakage while washing.
- Wide neck



Catalog No.	Capacity	Grad. Range	To Fit Rubber Stopper No	Pack Type	Pack Price	Case Type	Case Price
232.202.01	125 mL	50-125 mL	6	12	\$ 47.88	48	\$ 191.5
232.202.02	250 mL	50-200 mL	8	12	\$ 71.76	48	\$ 237.0
232.202.04	1000 mL	250-1000 mL	11	6	\$ 63.24	24	\$ 252.9

## Erlenmeyer Flask, Narrow Neck With Graduation

- Manufactured from ASTM E438 Type 1 Class A borosilicate 33 expansion glass.
- Feature durable white ceramic enamel scale which indicate approximate volume for measuring and mixing.
- Reinforced rim for easy pouring and also avoids any breakage while washing.



Catalog No.	Capacity	Grad. Range	To Fit Rubber Stopper No	Pack Type	Pack Price	Case Type	Case Price
231.202.01A	10 mL	5-10 mL	0	12	\$ 57.60	48	\$ 230.40
231.202.01	25 mL	10-20 mL	0	12	\$ 39.48	48	\$ 157.98
231.202.05	250 mL	50-200 mL	6	12	\$ 37.68	48	\$ 150.70
231.202.06	500 mL	100-500 mL	7	6	\$ 24.72	36	\$ 148.30
231.202.07	1000 mL	250-1000 mL	9	6	\$ 38.90	24	\$ 155.60
231.202.08	2000 mL	600-2000 mL	10	4	\$ 47.60	8	\$ 95.02
231.202.09	4000 mL	1000-4000 mL	10	Each	\$ 48.90	2	\$ 97.80
231.202.09A	5000 mL	1000-5000 mL	10	Each	\$ 85.90	2	\$ 171.80

**Sati International**

965 Reverchon St-Laurent (Quebec)  
Canada H4T 4L2

1 (833) 636-7850  
sales@satiinternational.com  
www.satiinternational.com