



WAXFILM

- TLG® Waxfilm are flexible, semi-transparent thermoplastic film with a two-way stretch
- Despite off-gassing, it creates a tight moisture barrier that maintains a strong seal and does not release from glassware.
- Non-toxic and free of plasticizers, primarily composed of polyolefins and paraffin waxes, with good resistance to acids, alkalis, and some organic solvents—use with acetone, diethyl ether, or tetrachloride is not recommended.
- Reduces the risk of tampered data with an odorless, colorless, and semi-transparent covering.
- The material's flexibility and self-sealing properties allow it to cling to irregular shapes and surfaces.
- Many polar substances, such as saline solutions, inorganic acids, and alkaline solutions, are resistant for up to 48 hours.

Catalog No.	Description	Diameter	Length (ft.)
CS-665WF	Waxfilm 2" x 250 ft.	2 inches	250
CS-666WF	Waxfilm 4 in.x125 ft.	4 inches	125