



Ready Stock  
Next Day Shipping

# SCINTILLATION VIALS

Limited Time  
**OFFER**  
Till Stock Last



These scintillation vials are produced from 33 expansion, low extractable Borosilicate glass conforming to USP Type I and ASTM E 438, Type I, Class A requirements. There are five trays per case and each tray contains 100 vials, shrink-wrapped to preserve cleanliness. The scintillation vials are supplied with unattached caps.

Catalog No.	Capacity	Height	Cap Size	Cap Type	Liner	Style	Case Size	Price
CS-1000-20	20.00 ml	57 mm	22-400	Polyethylene	Pulp-Backed Foil	Clear Glass	Case of 500	\$ 180.90/CS
CS-1200-20	20.00 ml	57 mm	24-400	Polyethylene	Pulp-Backed Foil	Clear Glass	Case of 500	\$ 180.90/CS

\*Minimum Order of 5 Cases of Each or Assorted.

\*All Prices are Fob Sati International Montreal.

 **VISIT US!** [WWW.SATIINTERNATIONAL.COM](http://WWW.SATIINTERNATIONAL.COM)

 **WRITE US!** [SALES@SATIINTERNATIONAL.COM](mailto:SALES@SATIINTERNATIONAL.COM)

 **CALL NOW**

**1 (833) 636-7850**