

TITRATION

SAMPLE BEAKERS



TLG® titration sample beakers are made from virgin polypropylene, have excellent visibility and very good resistance to acids.

- ▶ CS-MP-117 120 mL Titration Sample Beaker fits Metrohm USB and Robotic Sample Processors
- ▶ CS-MP-118 100 mL Titration Sample Beaker has graduations at 20, 40, 60, 80 and 100 mL and fits Mettler Toledo Autosamplers.
- ▶ CS-MP-129 120 mL sample beaker is a direct replacement for Metrohm 6.01400.200, Fits 602041030 sample rack for OMNIS Pick&Place

Catalog No.	Size/Description	Material	Bottom	Qty Per Case	Comparable to
CS-MP-117	40 x 113 mm (120 mL)	PP-CO	Flat	300	Metrohm™ 6.1459.300 & 6.1459.307
CS-MP-118	53.6 x 90.5mm (100 mL)	PP-H	Flat	450	Mettler Toledo 00101974, 511093388 & 51109389; Dyn-a-med 80099; Precision Weighing Balances 99008; Fisher Scientific NC9485544; VWR 56221-531 & 12578-108
CS-MP-129	44.7 x 113 mm (120 mL)	PP-CO	Flat	275	Metrohm™ 6.01400.200