



Worldwide scientific glassware manufacturer

ChemScience Centrifuge Tubes
(15/50 ML Tubes)

ChemScience Centrifuge Tubes (15/50 ML Tubes)

Where perfection is the prerequisite, even the smallest components get the same importance.

Even though centrifuge tubes are considered widely as just another part of the device, they are extremely important in the way they affect the outcome of the experiment.

Thus, at ChemScience, the highest standards are followed during the design & development of these tubes to ensure that there does not remain even the minutest chance of any flaws in this vital consumable.

As a result, these tubes of high quality are certified to be free from all affecting biological elements such as DNase, RNase, Endotoxin, Pyrogen and ATP.



15 ML



50 ML



**50 ML
Self Standing**

Perfection Is An Understatement

Experience the impact of perfection in your lab with ChemScience impeccable centrifuge tubes which are designed using the absolute best raw materials and procedures to ensure almost 100% purity. Moreover, these tubes are built to offer optimum convenience & durability.



Manufactured In Fully Automatic Human Touch-free Plant With Class 100k Clean Room Facility With Contolled Environment & Access



Produced From Swiss-made Diamond Finish Mold For Flawless 'Super Smooth' Internal Surface



Non Cytotoxic



Biologically Inert

Highly Transparent And Clear



US FDA Approved Virgin Polypropylene Without Any Trace Of Metals



Dnase, RNase, Endotoxin, Pyrogen Free



Can Withstand Temperature From -80 To +121°C



Perfection Is An Understatement



Easy Seal Design Of The Cap To Ensure Smooth Opening And Closing With Single Hand Repeatedly

Tamper Proof Graduation Marks & Writing Area - Resistant To Chemical & Physical Corrosion



Multiple Packaging Options

Fits All Standard And High Capacity Centrifuge Rotors



Clear Large Graduations On The Body Of The Tube



Superior Engineering Technique Enables Tubes To Withstand Higher G Force



15 ML
17000g



50 ML
20000g

Ordering Information:

ChemScience				
Free from Dnase/ Rnase/ Pyrogen/ Endotoxin				
Cat. No.	Description	Packing Type	Quantity / Pack	Quantity / Case
CS-CT-15-B	15 ML Centrifuge Tube With Flat Cap	Bulk	25 Tube	25 Tubes X 20 Packs
CS-CT-15-R	15 ML Centrifuge Tube With Flat Cap	Rack	25 Tube	25 Tubes X 10 Packs
CS-CT-15-B-S	15 ML Centrifuge Tube With Flat Cap, Sterile	Bulk	25 Tube	25 Tubes X 20 Packs
CS-CT-15-R-S	15 ML Centrifuge Tube With Flat Cap, Sterile	Rack	25 Tube	25 Tubes X 10 Packs
CS-CT-50-B	50 ML Centrifuge Tube With Flat Cap	Bulk	25 Tube	25 Tubes X 20 Packs
CS-CT-50-R	50 ML Centrifuge Tube With Flat Cap	Rack	25 Tube	25 Tubes X 20 Packs
CS-CT-50-B-S	50 ML Centrifuge Tube With Flat Cap, Sterile	Bulk	25 Tube	25 Tubes X 20 Packs
CS-CT-50-R-S	50 ML Centrifuge Tube With Flat Cap, Sterile	Rack	25 Tube	25 Tubes X 20 Packs
CS-CT-50-K-B	50 ML Self Standing Centrifuge Tube With Flat Cap	Bulk	25 Tube	25 Tubes X 20 Packs
CS-CT -50-K-B-S	50 ML Self Standing Centrifuge Tube With Flat Cap, Sterile	Bulk	25 Tube	25 Tubes X 20 Packs



ChemScience inc.
www.chemscience.com